



Intelligently Designed to Protect Augmented Bone




Geistlich Bio-Gide®

**Architect
the Future**

Your options when faced with these situations



Clinical case by Dr. Collin Campbell | United Kingdom



No Membrane?

Not using a membrane can lead to soft-tissue cell ingrowth into the graft material and slow vascularization with the consequence of more soft tissue and less bone volume.²¹⁻²⁴

Native Collagen Membrane!

Let's take a closer look on the role of a native collagen membrane

Cross-linked Collagen Membrane?

For certain defects, longer barrier functions are needed. Usually the collagen of these membranes are cross-linked.^{1,14} Several problems can be associated with cross-linked collagen membranes:



Causes foreign body reactions like recruitment of multinucleated giant cells^{1,2}



Causes inflammatory responses which leads to delayed healing, encapsulation², less tissue integration¹



Causes a high inflammation rate of soft-tissues within an early phase of healing (1 week) which is correlated with less bone volume after 6 months.⁴




A site treated with a cross-linked membrane shows a soft-tissue dehiscence 7 days post-surgery. Clinical case by Dr. Zumstein, Lucerne | Switzerland

Higher inflammation rate is linked to less bone volume for proper implant placement⁴

A Coordinated Interplay Between Vascularization and Barrier Time

Geistlich Bio-Gide® and Geistlich Bio-Oss® have demonstrated in numerous clinical studies sufficient gain of bone volume and quality.¹⁷⁻¹⁹




Setting the stage for optimal soft- and hard tissue regeneration^{1,12,16}

The most important requirements when choosing a suitable membrane for bone regeneration are biocompatibility, tissue integration, cell occlusiveness, nutrient transfer and handling.^{1,10,11} **Geistlich Bio-Gide® meets all these requirements.**

The natural collagen structure of Geistlich Bio-Gide® supports early vascularization resulting in good tissue integration and wound stabilization.^{1,12}

The smooth layer orientated towards the soft-tissue serves as a scaffold for the attachment of fibroblasts, while the barrier function prevents the ingrowth of soft-tissue into the newly forming bone underneath.^{1,7,12-14} The rough, open-pored bottom side serves as a framework for bone forming cells like osteoblasts^{14,15}, nutrients and growth factors and allows ingrowth of blood vessels.^{1,12,15} **These characteristics even make the handling of Geistlich Bio-Gide® easy and convenient. You have the choice between different product types and shapes – matching your specific needs.**



Geistlich Bio-Gide® acts as a temporary barrier instead of unnecessary blockade

Once the protective function of Geistlich Bio-Gide® has been fulfilled the membrane resorbs. The cells in the area to be regenerated are pre-determined for their desired function and will develop accordingly.^{1,8,9} The soft tissues heal without scarring and largely without further complications.⁵ **Thus, Geistlich Bio-Gide® has an appropriate barrier time.**^{1,5-9}



Be the architect of your patient's bone regeneration:

With Geistlich Bio-Gide® there is an ideally suited collagen membrane available to guide your daily regenerative needs^{1,24,10}



High long-term implant survival rate²⁰



After 12–14 years implants that were placed in bones augmented with Geistlich Bio-Gide® and Geistlich Bio-Oss® show the similar survival rate as implants in pristine bone²⁰



Reliable esthetic results: with using Geistlich Bio-Gide® and Geistlich Bio-Oss® for contour augmentation the facial bone wall was preserved in 95% of patients over a period of 5–9 years²⁰



More details about our distribution partners:
www.geistlich-biomaterials.com

Geistlich Pharma AG
Business Unit Biomaterials
Bahnhofstrasse 40
6110 Wolhusen, Schweiz
Phone +41 41 492 55 55
Fax +41 41 492 56 39
www.geistlich-biomaterials.com

Affiliate Australia and New Zealand
Geistlich Pharma Australia and New Zealand
The Zenith – Tower A,
Level 19, Suite 19.01
821 Pacific Highway
NSW 2067 Chatswood, Australia
Phone +61 1800 776 326
Fax +61 1800 709 698
info@geistlich.com.au
www.geistlich.com.au

Affiliate Great Britain and Ireland
Geistlich Sons Limited
1st Floor, Thorley House
Bailey Lane
Manchester Airport
Manchester M90 4AB, Great Britain
Phone +44 161 490 2038
Fax +44 161 498 6988
info@geistlich.co.uk
www.geistlich.co.uk


Affiliate North America
Geistlich Pharma North America Inc.
902 Carnegie Center, Suite 360
Princeton, NJ 08540 USA
Phone +1 855 799 5500
info@geistlich-na.com
www.geistlich-na.com

Distribution Canada
HANSamed Ltd.
2830 Argentia Road
Unit 5-8
L5N 8G4 Mississauga, Canada
Phone +1 800 363 2876
Fax +1 800 863 3213
orders@hansamed.net
www.hansamed.net

Even in such a difficult esthetic case including horizontal and vertical bone and tissue loss associated with recession on teeth 12 and 21, the combination of Geistlich Geistlich Bio-Oss® plus autologous bone and Geistlich Bio-Gide® provided a stable mucosal result around implant restoration for the long-term success.



Clinical case by Dr. Collin Campbell | United Kingdom



- 1 Initial situation: recession of tooth 12 and 21 are present
- 2 Retraction of flap shows considerable bone loss associated with tooth 12 and early bone loss associated with tooth 21. Vertical and horizontal defects associated with 11 implant site also clearly visible.
- 3 Additional picture with surgical guide in position demonstrating correct vertical position of implant.
- 4 Implant placed with cover screw. Autologous bone chips have been applied which were harvested locally.
- 5 Geistlich Bio-Oss® granules were applied over autologous bone to provide stability to the whole bone graft.
- 6 A double-layer of Geistlich Bio-Gide® was applied to protect the bone graft during healing period.
- 7 Immediate post-operative view following application of sutures.
- 8 2-year follow-up with a stable gingival position.



Geistlich Bio-Oss® plus autologous bone chips



Geistlich Bio-Gide®

References

- 1 Rothamel D et al. Clin. Oral Implants Res. 2005; 16(3): 369-378. (Pre-clinical study)
- 2 Al Maawi S et al. Semin Immunol. 2017 Feb;29:49-61. (Pre-clinical study)
- 3 Fickl S et al. Int J Periodontics Restorative Dent. 2018 Jan/Feb;38(1):e1-e7. (Clinical study)
- 4 Scheyer ET et al. J Clin Periodontol. 2016 Dec;43(12):1188-1199. (Clinical study)
- 5 Becker J et al. Clin Oral Implants Res. 2009; 20(7):742-749. (Clinical study)
- 6 Kim M et al. In Vivo. 2008; 22(2):231-6. (Pre-clinical study)
- 7 Zitzmann NU et al. Int J Oral Maxillofac Implants. 2012; 19(7):844-852. (Clinical study)
- 8 Giellkens PFM et al. Clin. Oral Implants Res. 2008;19:516-521. (Pre-clinical study)
- 9 Rothamel D et al. Int J Oral Maxillofac Implants. 2012 Jan–Feb;27(1):146–54. (Pre-clinical study)
- 10 Caballé-Serrano J et al. J Clin Exp Dent. 2018 May 1;10(5):e477-e483. (Systematic review)
- 11 Benic GI, Hämmerle CH. Periodontol 2000. 2014 Oct;66(1):13-40. (Systematic review)
- 12 Schwarz F et al. Clin. Oral Implants Res. 2006;17(4):403-409. (Pre-clinical study)
- 13 Tal H et al. Clin Oral Implants Res. 2008; 19(3) : 295-302. (Clinical study)
- 14 Rothamel D et al. Clin. Oral Implants Res. 2004;15:443-449. (Pre-clinical study)
- 15 Data on File. Geistlich Pharma AG, Wolhusen, Switzerland. (Non-clinical)
- 16 Schwarz F et al. Clin Oral Implants Res. 2014 Sept;25(9):1010-5. (Clinical study)
- 17 Wallace SS et al. Int J Periodontics Restorative Dent. 2005;25(6):551-559. (Clinical study)
- 18 Wallace SS et al. Ann Periodontol. 2003;8(1):328-43. (Clinical study)
- 19 Norton MR. Clin Oral Implants Res. 2001 Feb;12(1):79-84. (Clinical study)
- 20 Jung RE et al. Clin. Oral Implants Res. 2013 Oct;24(10):1065-73. (Clinical study)
- 21 Dahlin C. Plast Reconstr Surg. 1988 May;81(5):672-6. (Clinical study)
- 22 Omar O et al. J Clin Periodontol. 2019 Jun;46 Suppl 21(Suppl Suppl 21):103-123. (Review)
- 23 Schenk RK et al. Int J Oral Maxillofac Implants. Jan-Feb 1994;9(1):13-29. (Pre-clinical study)
- 24 Perelman-Karmon M et al. Int J Periodontics Restorative Dent. 2012 Aug;32(4):459-65. (Clinical study)