Extraction Socket Management

Indication Sheet | Prof.Ki-Tae Koo



Indication Sheet
Extraction Socket Management

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Treatment Concepts for Infected Extraction Sockets

by Pf. Ki-Tae Koo, Seoul National University,
Department of Periodontology, Seoul, Korea

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
Preservation & membrane application

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

Classification of Infected Extraction Sockets

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Surgical Procedure

Type I

Defect configuration	4-wall socket / No soft tissue loss
Healing potential	Good
Pathology	Endodontic origin, fracture, severe caries
Type of bone graft	Particulated or matrix type
Soft tissue seal	흡수성 콜라겐 멤브레인 or 연조직 대체제
Timing of implant placement	Immediate (Ridge preservation 시 delayed 시행)
Predictability	Highly predictable

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

Type II

Defect configuration	3-wall socket / No soft tissue loss
Healing potential	Good
Pathology	Periodontal origin, Periodontal-Endodontic combined
Type of bone graft	Matrix type
Soft tissue seal	흡수성 콜라겐 멤브레인
Timing of implant placement	Immediate or early
Predictability	Slightly compromised

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



Type III

Defect configuration	3-walled socket / Soft tissue loss
Healing potential	Poor
Pathology	Periodontal origin, Periodontal-Endodontic combined
Type of bone graft	Matrix type
Soft tissue seal	흡수성 콜라겐 멤브레인
Timing of implant placement	Delayed
Predictability	Highly compromised soft tissue augmentation necessary

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
Type IV-A


 <p>Type IV-A</p>	Defect configuration	4-wall or 2-wall socket / No soft tissue loss
	Healing potential	Poor
	Pathology	Periodontal origin, Periodontal-Endodontic combined
	Type of bone graft	Matrix type
	Soft tissue seal	흡수성 콜라겐 멤브레인
	Timing of implant placement	Delayed
	Predictability	Highly compromised

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

Type V


 <p>Type V</p>	Defect configuration	4-wall or 2-wall socket / soft tissue loss
	Healing potential	Poor
	Pathology	Periodontal origin, Periodontal-Endodontic combined
	Type of bone graft	Particulated or Matrix type
	Soft tissue seal	흡수성 콜라겐 멤브레인
	Timing of implant placement	Delayed
	Predictability	Compromised vertical augmentation necessary

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

Type IV-B

 <p>Type IV-B</p>	Defect configuration	2-wall socket / No soft tissue loss
	Healing potential	Very poor
	Pathology	Periodontal origin, Periodontal-Endodontic combined
	Type of bone graft	Matrix type
	Soft tissue seal	흡수성 콜라겐 멤브레인
	Timing of implant placement	Delayed
	Predictability	Highly compromised

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
Geistlich
Bio-Gide®
천연 콜라겐 차폐막

LITERATURE REFERENCES

- Ahn, J. J. & Shin, H. I. (2008) Bone tissue formation in extraction sockets from sites with advanced periodontal disease: A histomorphometric study in humans. *International Journal of Oral Maxillofacial Implants* 23: 1133-1138.
- Araujo, M., Linder, E., Wennstrom, J. & Lindhe, J. (2008) The influence of bio-oss collagen on healing of an extraction socket: An experimental study in the dog. *International Journal of Periodontics & Restorative Dentistry* 28:123-135.
- Araujo, M. G., Carmagnola, D., Berglundh, T., Thilander, B. & Lindhe, J. (2001) Orthodontic movement in bone defects augmented with bio-oss. An experimental study in dogs. *Journal of Clinical Periodontology* 28: 73-80.
- Araujo, M. G., Liljenberg, B. & Lindhe, J. (2010) Dynamics of bio-oss collagen incorporation in fresh extraction wounds: An experimental study in the dog. *Clinical Oral Implants Research* 21: 55-64.
- Araujo, M. G. & Lindhe, J. (2005) Dimensional ridge alterations following tooth extraction. An experimental study in the dog. *Journal of Clinical Periodontology* 32: 212-218.
- Caplanis, N., Lozada, J. L. & Kan, J. (2009) Extraction defect: Assessment, classification and management. *International Journal of Clinical Implant Dentistry* 1:1-11.
- Cardaropoli, D., Tamagnone, L., Roffredo, A., Gaviglio, L. & Cardaropoli, G. (2012) Socket preservation using bovine bone mineral and collagen membrane: A randomized controlled clinical trial with histologic analysis. *International Journal of Periodontics & Restorative Dentistry* 32: 421-430.
- Cardaropoli, G., Araujo, M., Hayacibara, R., Sukekava, F. & Lindhe, J. (2005) Healing of extraction sockets and surgically produced - augmented and non-augmented - defects in the alveolar ridge. An experimental study in the dog. *Journal of Clinical Periodontology* 32: 435-440.
- Cardaropoli, G., Araujo, M. & Lindhe, J. (2003) Dynamics of bone tissue formation in tooth extraction sites. An experimental study in dogs. *Journal of Clinical Periodontology* 30: 809-818.
- Davies, J. E. (1998) Mechanisms of endosseous integration. *The International Journal of Prosthodontics* 11: 391- 401.
- Davies, J. E. (2003) Understanding peri-implant endosseous healing. *Journal of Dental Education* 67: 932-949.
- Fanuscu, M. I. & Chang, T. L. (2004) Three-dimensional morphometric analysis of human cadaver bone: microstructural data from maxilla and mandible. *Clinical Oral Implants Research* 15: 213-218.
- Elian, N., Cho, S. C., Froum, S., Smith, R. & Tarnow, D. (2007) A simplified socket classification and repair technique. *Practical Procedures and Aesthetic Dentistry* 19:99-104.
- Hahn, M., Vogel, M., Pompesius-Kempa, M. & Delling, G. (1992) Trabecular bone pattern factor--a new parameter for simple quantification of bone microarchitecture. *Bone* 13: 327-330.
- Hammerle, C. H., Chiantella, G. C., Karring, T. & Lang, N. P. (1998) The effect of a deproteinized bovine bone mineral on bone regeneration around titanium dental implants. *Clinical Oral Implants Research* 9: 151-162.
- Hildebrand, T. & Rueggsegger, P. (1997) Quantification of bone microarchitecture with the structure model index. *Computer Methods in Biomechanics and Biomedical Engineering* 1: 15-23.
- Kilkenny, C., Browne, W., Cuthill, I. C., Emerson, M., Altman, D. G., National Centre for the Replacement, R. & Reduction of Animals in, R. (2011) Animal research: Reporting in vivo experiments--the arrive guidelines. *Journal of Cerebral Blood Flow & Metabolism* 31: 991-993.
- Kim, J. H., Susin, C., Min, J. H., Suh, H. Y., Sang, E. J., Ku, Y., Wikesjo, U. M. & Koo, K. T. (2014) Extraction sockets: Erratic healing impeding factors. *Journal of Clinical Periodontology* 41: 80-85.
- Lindhe, J., Cecchinato, D., Donati, M., Tomasi, C. & Liljenberg, B. (2014) Ridge preservation with the use of deproteinized bovine bone mineral. *Clinical Oral Implants Research* 25: 786-790.
- Lopez-Martinez, F., Gomez Moreno, G., Olivares-Ponce, P., Eduardo Jaramillo, D., Eduardo Mate Sanchez de Val, J. & Calvo-Guirado, J. L. (2015) Implants failures related to endodontic treatment. An observational retrospective study. *Clinical Oral Implants Research* 26: 992-995.
- Marcaccini, A. M., Novaes, A. B., Jr., Souza, S. L., Taba, M., Jr. & Grisi, M. F. (2003a) Immediate placement of implants into periodontally infected sites in dogs. Part 2: A fluorescence microscopy study. *International Journal of Oral Maxillofacial Implants* 18: 812-819.
- Moon, H. S., Won, Y. Y., Kim, K. D., Ruprecht, A., Kim, H. J., Kook, H. K. & Chung, M. K. (2004) The three dimensional microstructure of the trabecular bone in the mandible. *Surgical and Radiologic Anatomy* 26: 466-473.
- Parfitt, A. M., Drezner, M. K., Glorieux, F. H., Kanis, J. A., Malluche, H., Meunier, P. J., Ott, S. M. & Recker, R. R. (1987) Bone histomorphometry: standardization of nomenclature, symbols, and units. Report of the ASBMR Histomorphometry Nomenclature Committee. *Journal of Bone and Mineral Research* 2: 595-610.
- Pietrokovski, J. & Massler, M. (1967) Alveolar ridge resorption following tooth extraction. *Journal of Prosthetic Dentistry* 17: 21-27.
- Roze, J., Babu, S., Saffarzadeh, A., Gayet-Delacroix, M., Hoornaert, A. & Layrolle, P. (2009) Correlating implant stability to bone structure. *Clinical Oral Implants Research* 20: 1140-1145.
- Schroeder, H. E. & Munzel-Pedrazzoli, S. (1973) Correlated morphometric and biochemical analysis of gingival tissue. Morphometric model, tissue sampling and test of stereologic procedures. *Journal of microscopy* 99: 301-329.
- Trombelli, L., Farina, R., Marzola, A., Bozzi, L., Liljenberg, B. & Lindhe, J. (2008) Modeling and remodeling of human extraction sockets. *Journal of Clinical Periodontology* 35: 630-639.

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